CONTRAINDICATIONS

There are no contraindications for the use of AJL BLUE when used according to these instructions.

PRECAUTIONS

It is recommended that after injection any excess AJL BLUE is removed immediately from the eye by irrigation/aspiration techniques with saline solution.

Known risks of a possible reuse of the product are:

- Reusing this product creates a potential risk of infection for the patient. Contamination of the product may cause injury to the patient.
- Reprocessing of the product may compromise its structural integrity, as well as its sterility and functionality.
- It should not be used in patients with hypersensitivity.
- It is not recommended for use in children and pregnant women.

ADVERSE REACTIONS

Inadvertent staining of the posterior capsule of the lens and the vitreous face. Staining of the vitreous capsule is usually self-limiting, lasting up to one week.

May cause allergic reactions in hypersensitive patients.

There may be cases of postoperative inflammatory reactions and corneal edema.

CLINICAL APPLICATIONS

AJL BLUE selectively stains connective tissue structures in the human eye, such as the anterior lens capsule.

AJL BLUE is designed to be applied directly in the anterior lens capsule, staining any part of the capsule that comes into contact with the dye. The excess dye is removed by aspiration and/or irrigation with saline solution.

The dye does not penetrate the posterior capsule, permitting visualization of the anterior capsule by contrast with the lens and unstained cortex thereof.

INSTRUCTIONS FOR USE

Open the pack and remove the contents within a sterile operating field. Connect the luer lock cannula to the end of the syringe and confirm correct assembly. Press down on the plunger gently to expel a few drops of the product to prevent the introduction of air bubbles into the anterior chamber. The syringe is now ready for use.

SURGICAL PROCEDURE

The surgical procedure is the responsibility of each surgeon, who must be previously trained in the use of this surgical technique. The use of the 27G cannula included with the product is recommended, as is, the removal of excess dye by means of aspiration and/or irrigation techniques with saline solution.

The manufacturer is not liable for any damage caused to the patient resulting from the surgical technique used by the surgeon to apply the product.

Reporting customer complaints including quality complaints, adverse events and other medical device related observations:







Customer complaints including quality complaints, adverse events and other medical device related observations should be reported to AJL without delay. A report describing the details of the complaint/event, the applied therapy, the product type, LOT/serial number of the medical device used is requested.

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